



# Beyond the Incubator: A Comprehensive Approach to Contamination Control for Specialty Equipment and Connected Environments

*Contamination Control Strategies for Complex Environments*



Elizabeth Massoth  
*Biocontainment Specialist*  
CURIS System

## Learning Objectives

1. **Assess whether contamination originating in specialty equipment is localized or has spread to connected cleanroom environments and HVAC systems.**
2. **Compare decontamination approaches**—including manual, OEM, and gaseous biodecontamination technologies—and evaluate their limitations for complex equipment geometries.
3. **Apply a tiered contamination response model** (equipment → room → HVAC) appropriate to scope, contamination vector, and regulatory expectations.
4. **Mitigate risk** across mold, bacterial, viral, and DNA fragmentation contamination vectors in specialty equipment environments.

## Why This Matters: Regulatory & Business Stakes

### Regulatory Stakes

- **Annex 1 (2022):** equipment is explicitly within CCS scope
- **Inspectors expect documented** equipment-level contamination risk assessments
- **ICH Q9:** risk-based justification required for response decisions
- **483 / Warning Letter exposure** when non-traditional sources aren't addressed in the CCS

### Operational & Business Stakes


- **Line downtime**—lost batches every day operations are paused
- **Deviation & CAPA cycle time**—investigation, root cause, remediation
- **Sterility assurance & batch disposition** consequences
- **OEM service scheduling** can extend downtime by weeks
- **Cross-product / campaign contamination risk** in multi-product facilities

## When Equipment Becomes the Source: The Contamination Cascade

**Contamination rarely stays contained to a single piece of equipment**

- Airflow Dynamics
- Human Cross-Contamination
- Equipment Design
- HVAC Integration

**The result: an isolated equipment event becomes a facility-wide incident requiring costly and disruptive intervention**



## Real-World Scenarios: When Contamination Escalates



### Incubators

Mold migrating to adjacent cleanroom surfaces via door opening and airflow



### Cell Sorters

Aerosol-generating operation creates room-level contamination, triggering unplanned full-room decontamination



### Lyophilizers

Internal bioburden in chambers and condenser manifolds compromising batch sterility assurance



### Depyrogenation Tunnels

High-temperature environment and inaccessible internal surfaces create persistent contamination risk

### Common Pattern Across Events:

**late detection → inadequate initial response scope → facility-wide intervention required**

# Decontamination Methods

*Evaluating Your Options for Specialty Equipment*

## Why “Clean the Equipment and Move On” Often Falls Short

*Three options dominate the field: manual spray-and-wipe, OEM service, and gaseous/vaporous biodecontamination. The first two are the most common—and the most likely to leave contamination behind.*




**Manual surface methods can't reach the inside** Spray-and-wipe addresses external surfaces only. Technicians can't reliably reach internal components, HEPA housings, condensate pathways, or hard-to-access areas—and efficacy is highly operator-dependent.

**OEM service is partial by design** Service interventions are equipment-specific and typically address the equipment only—not the surrounding room or shared airflow. Scheduling constraints can extend downtime by weeks.




**Disinfectant choice carries its own risks** Insufficient kill claims (lack of EPA testing/approval), misapplied contact times, and chemistries that leave residues—harmful to staff, products, or capable of creating germ reservoirs. Repeat contamination events typically trace back to one of these gaps.

## Spaulding Hierarchy & Log Reduction: Understanding Disinfection Difficulty

<p>More Resistant</p> <p>Less Resistant</p>	<b>Bacterial Spores</b>	Spores of: <i>C. difficile</i> ; <i>C. tetanus</i> ; <i>Geobacillus Stearothermophilus</i> ; <i>C. botulinum</i> ; <i>C. sporogenes</i> ; <i>Bacillus anthracis</i> ; <i>Bacillus atrophaeus</i>
	<b>Mycobacteria</b>	<i>M. tuberculosis</i> (TB); <i>M. avian</i>
	<b>Viruses Without Envelopes</b>	Norovirus; Rotavirus; Rhinovirus; Poliovirus; Papillomavirus (HPV); Adenovirus x 10 <sup>10</sup> ; Coxsackievirus; Calicivirus
	<b>Fungi including Fungal Spores</b>	<i>Aspergillus fumigatus</i> ; <i>A. flavus</i> ; <i>A. niger</i> ; <i>Candida albicans</i> ; <i>Calonectria henricotiae</i> ; <i>Calonectria pseudonaviculata</i>
	<b>Gram Negative Bacteria</b>	<i>Pseudomonas</i> ; <i>Acinetobacter</i> ; <i>E. coli</i> ; Enterobacteriaceae; <i>Legionella</i> ; <i>Salmonella</i>
	<b>Gram Positive Bacteria</b>	<i>Staphylococcus</i> ; <i>Enterococcus</i> (VRE); <i>Streptococcus</i> ; Clostridia vegetative rods; <i>Rathayibacter agropyri</i> ; <i>Rathayibacter iranicus</i> , <i>C. Bovis</i>
	<b>Viruses With Lipid Envelopes</b>	Influenza; HBV; HCV; HIV; RSV; Coronavirus; CMV; HSV; Measles; Mumps; Rubella; VZV (Varicella-Zoster) Shingles/Chickenpox



## Evaluating Chemistry & Form Factor

		Chlorine Dioxide	35% Hydrogen Peroxide	Ionized Hydrogen Peroxide	7% Hybrid Hydrogen Peroxide™
Material Compatibility		●	●	●	●
Staff & Animal Safety (exposure to leakage / residues)		●	●	●	●
Residue Profile		●	●	●	●
Portability & Setup		●	●	●	●
Multi-Tier Flexibility		●	●	●	●
● Favorable		● Conditional / Limited		● Concern	
Selection should be driven by contamination vector, equipment materials, and required contact time Verify EPA registration, sporicidal claims, and material compatibility data for your specific application.					

## Internal Geometries & Material Compatibility



### Particle Counters

Laser optics in sensor housing—sensitive to oxidizing agents and moisture.



### Lyophilizers

Internal shelving, manifolds, condenser coils—seals must withstand gaseous decon conditions.



### Cell Sorters

Fluidics paths, sort blocks, optics—complex geometry; gaskets and tubing require compatibility assessment.



### Depyrogenation Tunnels

Inaccessible high-temperature internals—evaluate conveyor and seal integrity.

*Across all four:*

- Vapor distribution within sealed enclosures without direct injection
- BIs and CIs at worst-case internal positions to confirm coverage
- Cycle qualification specific to equipment type and configuration
- Lower-concentration chemistries preserve sensitive electronics and optics


*CURIS field research:* HHP™ compatibility validated on particle counters and cleanroom surfaces.

# A Tiered Approach to Contamination Control

*Equipment → Room + Equipment → HVAC & Surrounding Spaces / Whole Facility*

## The Tiered Model: Scope Determines Strategy

1. Address the affected equipment directly.
  2. Expand to the room, adjacent spaces (MALs/PALs), and shared airflow paths
  3. Full-scale response involving multiple zones, HVAC systems, and potentially professional services.
- Scale of response should be driven by EM data, contamination type, and facility risk assessment—not assumption
  - CCS documentation should define the trigger criteria and decision logic for each level of response
  - Manual disinfection is labor-intensive, imperfect, and hard to validate. Consider where easy-to-implement vapor decon could measurably improve your contamination control efforts.



## Regulatory Alignment: Annex 1, cGMP, and CCS Expectations

- **EU GMP Annex 1 (2022): CCS must address all sources of contamination—equipment is explicitly within scope**
- Risk-based approach required: document your assessment, decision criteria, technology rationale, and remediation outcomes
- **The tiered model supports inspection readiness:**
  - Clear definition of scope at each response level
  - Validated decontamination approach with documented efficacy data
  - EM data as the primary validation outcome for remediation effectiveness
  - Ongoing monitoring program demonstrates sustained contamination control
- **Treating contamination as a facility-wide ecosystem problem—not just an isolated equipment issue—aligns with modern CCS philosophy**

---

# Field Validation: Specialized Equipment Fumigation

## Decontaminating Cell Sorters in ATMP Facilities

- **Setting**




Cell sorters can require decon at maintenance, calibration, and campaign changeover.

- **Challenge**

Complex geometries. **Sensitive optics and electronics.**

- **Solution**

Hybrid hydrogen peroxide, sorter in place—no disassembly.  **$\geq 6$ -log kill confirmed at worst-case positions.** Return to service without OEM service call.



## Beyond Microbial Kill: DNA Amplicon Decontamination

- **Setting**

PCR-based testing in CGT and QC environments. **Trace DNA amplicons cause false positives.**

- **Challenge**

Thermal cyclers and PCR instruments can't tolerate high-consequence / high-concentration DNA-degrading chemistries. Amplicons on internal airflow components are hard to reach. **Replacement is the common (COSTLY) fallback.**

- **Solution**

7% hybrid hydrogen peroxide, pulsed delivery, sealed BSC. Pre-treatment: detectable amplicons. **Post-treatment: zero signal at 18/18 locations for all 3 test cycles.** Equipment integrity preserved.



Photo 1




Photo 10

## Treating What You Can't Reach: Interstitial Space and Filter Banks

### Interstitial Space

**Setting** Above-ceiling plenum implicated as mold and particulate reservoir. Inaccessible to manual cleaning.

**Challenge** Staff can't reach interstitial space. Manual disinfection not possible.


**Solution** Room and plenum fumigated simultaneously — ceiling tiles removed, fans directed upward through apertures. Validated kill at worst-case BI positions throughout both volumes.

### HEPA Filter Banks

**Setting** HEPA filter banks implicated by EM data. Changeout would mean an extended shutdown and requalification.

**Challenge** Filters can't be wiped or visually verified. Need validated kill—without changeout.

**Solution** Gaseous decon delivered through HVAC.  $\geq 6$ -log kill confirmed at upstream, in-housing, and downstream BIs.



Same principle: validated gaseous kill where manual cleaning can't reach.

## Bulk Autoclave Repurposed for Material Fumigation

### • Setting


Getinge Bulk Autoclave for **material transfer**. Heat- and moisture-sensitive items also required a validated transfer pathway.

### • Challenge

**Sensitive items cannot tolerate autoclaving.** No dedicated room available for a separate material fumigation chamber.

### • Solution

Existing autoclave **converted to dual-purpose fumigation chamber** via decouplable hybrid hydrogen peroxide system with onboard drying. Autoclave retained full sterilization functionality. Validated kill at worst-case BI positions.



## Key Takeaways

- **Specialty equipment is a significant and underappreciated contamination source—incubators, cell sorters, lyophilizers, and similar devices can serve as reservoirs that seed broader facility contamination**
- **Contamination rarely stays isolated: airflow, human handling, and HVAC integration create pathways that must be systematically addressed**
- **No single method covers all scenarios: technology selection should be driven by contamination vector, equipment geometry, material compatibility, and required validation rigor**
- **The tiered approach—equipment → room → HVAC—provides a scalable, evidence-based framework for contamination response**
- **Align contamination response with your CCS: documented scope, validated outcomes, and EM data are the foundation of regulatory readiness**

# Thank You

*Questions & Discussion*

*[beth.massoth@curissystem.com](mailto:beth.massoth@curissystem.com)*